



The new FLOM GASTORR BG series

- Innovative
- new vacuum pump
- variable vacuum
- easy maintenance



Product name	Gastorr			
Part number	BG-32-02	BG-34-02	BG-32-10	BG-34-10
Application	Analytical		Semi-prep	
Flow rate	2mL/min		10mL/min	
Vacuum Level	70~170hPa (adjustable)			
Channels	2	4	2	4
Internal volume (1 channel)	510uL		1021uL	
Wetted materials	Teflon AF®, ETFE, PPS			
Communication	RS232C			
(WxDxHmm)	70 X 270 X 130			
Power consumption	0.8A, 80VA			
Voltage range	AC100~240V 50/60Hz			
Weight	2.4kg			

