



The new FLOM INERT PUMP series

- Innovative
- very stable flow
- chemical inert
- high dosing pressure



Product name	Dual Pump, inert	
Part number	KP-22-33D	KP-22-13D
Application	Semi-prep	Analytical
Max pressure	2 MPa	2 MPa
Flow rate	0.1 – 50 mL/min	0.01 – 10.0 mL/min
1 head volume	154 uL	32 uL
Flow stability	+/-0.3%	+/-0.3%
Features	Pressure pulse reduction, flow volume calibration, pressure zero reset, dispensing mode, upper pressure limiter	Pressure pulse reduction, flow volume calibration, pressure zero reset, dispensing mode, upper pressure limiter
Pressure accuracy	+/-5%	+/-5%
Wetted materials	PTFE, ETFE, Chemraz	PTFE, ETFE, Chemraz
Remote input	Pump on/off	Pump on/off
Remote output	Pump error	Pump error
Serial port	RS232C	RS232C
Dimensions	144H x 79W x 285D	144H x 79W x 285D
Weight	3.8 kg	3.8 kg

